

SETD3

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 67kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SETD3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C14orf154;

To place an order, please [Click HERE](#).

Catalog #: A8071

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84193

Isotype: IgG

Swiss Prot: Q86TU7

Purity: Affinity purification

For research use only.