

LHX4

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LHX4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CPHD4;

Catalog #: A8072

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 89884

Isotype: IgG

Swiss Prot: Q969G2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a large protein family which contains the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein is a transcription factor involved in the control of differentiation and development of the pituitary gland. Mutations in this gene cause combined pituitary hormone deficiency 4.

To place an order, please [Click HERE](#).