

ADAMTSL1

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 193kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ADAMTSL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1;

Catalog #: A8073

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 92949

Isotype: IgG

Swiss Prot: Q8N6G6

Purity: Affinity purification

For research use only.

Background:

This gene encodes a secreted protein and member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family. This protein lacks the metalloproteinase and disintegrin-like domains, which are typical of the ADAMTS family, but contains other ADAMTS domains, including the thrombospondin type 1 motif. This protein may have important functions in the extracellular matrix. Alternative splicing results in multiple transcript variants encoding distinct proteins.

To place an order, please [Click HERE](#).