

MonoMethyl-NONO-R357/R365

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human NONO

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

P54; NMT55; NRB54; P54NRB; PPP1R114;

Catalog #: A8087

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4841

Isotype: IgG

Swiss Prot: Q15233

Purity: Affinity purification

For research use only.

Background:

This gene encodes an RNA-binding protein which plays various roles in the nucleus, including transcriptional regulation and RNA splicing. A rearrangement between this gene and the transcription factor E3 gene has been observed in papillary renal cell carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes exist on Chromosomes 2 and 16.

To place an order, please [Click HERE](#).