

IL1RL1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 63kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL1RL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

T1; ST2; DER4; ST2L; ST2V; FIT-1; IL33R;

Catalog #: A8097

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9173

Isotype: IgG

Swiss Prot: Q01638

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the interleukin 1 receptor family. Studies of the similar gene in mouse suggested that this receptor can be induced by proinflammatory stimuli, and may be involved in the function of helper T cells. This gene, interleukin 1 receptor, type I (IL1R1), interleukin 1 receptor, type II (IL1R2) and interleukin 1 receptor-like 2 (IL1RL2) form a cytokine receptor gene cluster in a region mapped to chromosome 2q12. Alternative splicing of this gene results in multiple transcript variants.

To place an order, please [Click HERE](#).