

## ACRV1

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**28kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human ACRV1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

SP-10; SPACA2; D11S4365;

**Catalog #:**A8098

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**56

**Isotype:**IgG

**Swiss Prot:**P26436

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a testis-specific, differentiation antigen, acrosomal vesicle protein 1, that arises within the acrosomal vesicle during spermatogenesis, and is associated with the acrosomal membranes and matrix of mature sperm. The acrosomal vesicle protein 1 may be involved in sperm-zona binding or penetration. Alternatively spliced transcript variants have been described.

*To place an order, please [Click HERE](#).*