

ARSF

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:66kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ARSF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ASF;

Catalog #:A8102

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:416

Isotype:IgG

Swiss Prot:P54793

Purity:Affinity purification

For research use only.

Background:

This gene is a member of the sulfatase family, and more specifically, the arylsulfatase subfamily. Members of the subfamily share similarity in sequence and splice sites, and are clustered together on chromosome X, suggesting that they are derived from recent gene duplication events. Sulfatases are essential for the correct composition of bone and cartilage matrix. The activity of this protein, unlike that of arylsulfatase E, is not inhibited by warfarin. Multiple alternatively spliced variants, encoding the same protein, have been identified.

To place an order, please [Click HERE](#).