

ATP5O

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 23kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ATP5O

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ATPO; OSCP; HMC08D05;

Catalog #: A8103

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 539

Isotype: IgG

Swiss Prot: P48047

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a component of the F-type ATPase found in the mitochondrial matrix. F-type ATPases are composed of a catalytic core and a membrane proton channel. The encoded protein appears to be part of the connector linking these two components and may be involved in transmission of conformational changes or proton conductance.

To place an order, please [Click HERE](#).