

CD6

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CD6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TP120;

Catalog #: A8107

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 923

Isotype: IgG

Swiss Prot: P30203

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein found on the outer membrane of T-lymphocytes as well as some other immune cells. The encoded protein contains three scavenger receptor cysteine-rich (SRCR) domains and a binding site for an activated leukocyte cell adhesion molecule. The gene product is important for continuation of T cell activation. This gene may be associated with susceptibility to multiple sclerosis (PMID: 19525953, 21849685). Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).