

## CDH6

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 88kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human CDH6

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CAD6; KCAD;

**Catalog #:** A8109

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 1004

**Isotype:** IgG

**Swiss Prot:** P55285

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. The encoded protein is a type II cadherin and may play a role in kidney development as well as endometrium and placenta formation. Decreased expression of this gene may be associated with tumor growth and metastasis.

*To place an order, please [Click HERE](#).*