

## CSTF2

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:2000

**Calculated MW:** 61kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human CSTF2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CstF-64;

**Catalog #:** A8116

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 1478

**Isotype:** IgG

**Swiss Prot:** P33240

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a nuclear protein with an RRM (RNA recognition motif) domain. The protein is a member of the cleavage stimulation factor (CSTF) complex that is involved in the 3' end cleavage and polyadenylation of pre-mRNAs. Specifically, this protein binds GU-rich elements within the 3'-untranslated region of mRNAs.

*To place an order, please [Click HERE](#).*