

CFD

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human CFD

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DF; ADN; PFD; ADIPSIN;

Catalog #:A8117

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1675

Isotype:IgG

Swiss Prot:P00746

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the trypsin family of peptidases. The encoded protein is a component of the alternative complement pathway best known for its role in humoral suppression of infectious agents. This protein is also a serine protease that is secreted by adipocytes into the bloodstream. Finally, the encoded protein has a high level of expression in fat, suggesting a role for adipose tissue in immune system biology.

To place an order, please [Click HERE](#).