

LAMB3

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 130kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human LAMB3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AI1A; LAM5; LAMNB1; BM600-125KDA;

Catalog #: A8127

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3914

Isotype: IgG

Swiss Prot: Q13751

Purity: Affinity purification

For research use only.

Background:

The product encoded by this gene is a laminin that belongs to a family of basement membrane proteins. This protein is a beta subunit laminin, which together with an alpha and a gamma subunit, forms laminin-5. Mutations in this gene cause epidermolysis bullosa junctional Herlitz type, and generalized atrophic benign epidermolysis bullosa, diseases that are characterized by blistering of the skin. Multiple alternatively spliced transcript variants that encode the same protein have been found for this gene.

To place an order, please [Click HERE](#).