

NPR3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NPR3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NPRC; ANP-C; ANPRC; NPR-C; ANPR-C; GUCY2B; C5orf23;

Catalog #: A8138

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4883

Isotype: IgG

Swiss Prot: P17342

Purity: Affinity purification

For research use only.

Background:

This gene encodes one of three natriuretic peptide receptors. Natriuretic peptides are small peptides which regulate blood volume and pressure, pulmonary hypertension, and cardiac function as well as some metabolic and growth processes. The product of this gene encodes a natriuretic peptide receptor responsible for clearing circulating and extracellular natriuretic peptides through endocytosis of the receptor. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).