

STXBP3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 68kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STXBP3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PSP; MUNC18C; UNC-18C; MUNC18-3;

To place an order, please [Click HERE](#).

Catalog #: A8153

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6814

Isotype: IgG

Swiss Prot: O00186

Purity: Affinity purification

For research use only.