

UCHL3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 26kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UCHL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UCH-L3;

Catalog #: A8156

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7347

Isotype: IgG

Swiss Prot: P15374

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the deubiquitinating enzyme family. Members of this family are proteases that catalyze the removal of ubiquitin from polypeptides and are divided into five classes, depending on the mechanism of catalysis. This protein may hydrolyze the ubiquitinyl-N-epsilon amide bond of ubiquitinated proteins to regenerate ubiquitin for another catalytic cycle. Alternative splicing results in multiple transcript variants that encode different protein isoforms.

To place an order, please [Click HERE](#).