

SLC25A11

Reactivity: Human

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:100

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC25A11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OGC; SLC20A4;

Catalog #: A8163

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8402

Isotype: IgG

Swiss Prot: Q02978

Purity: Affinity purification

For research use only.

Background:

The oxoglutarate/malate carrier transports 2-oxoglutarate across the inner membranes of mitochondria in an electroneutral exchange for malate or other dicarboxylic acids (summary by Iacobazzi et al., 1992)

To place an order, please [Click HERE](#).