

STK16

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STK16

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KRCT; MPSK; TSF1; PKL12;

To place an order, please [Click HERE](#).

Catalog #: A8166

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8576

Isotype: IgG

Swiss Prot: O75716

Purity: Affinity purification

For research use only.