

STX11

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STX11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FHL4; HLH4; HPLH4;

Catalog #: A8169

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8676

Isotype: IgG

Swiss Prot: O75558

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the syntaxin family. Syntaxins have been implicated in the targeting and fusion of intracellular transport vesicles. This family member may regulate protein transport among late endosomes and the trans-Golgi network. Mutations in this gene have been associated with familial hemophagocytic lymphohistiocytosis.

To place an order, please [Click HERE](#).