

GPRC5A

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GPRC5A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAI3; TIG1; RAIG1; GPCR5A; PEIG-1;

Catalog #: A8173

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9052

Isotype: IgG

Swiss Prot: Q8NFJ5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the type 3 G protein-coupling receptor family, characterized by the signature 7-transmembrane domain motif. The encoded protein may be involved in interaction between retinoid acid and G protein signalling pathways. Retinoic acid plays a critical role in development, cellular growth, and differentiation. This gene may play a role in embryonic development and epithelial cell differentiation.

To place an order, please [Click HERE](#).