

THOC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human THOC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P84; HPR1; P84N5;

Catalog #: A8179

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9984

Isotype: IgG

Swiss Prot: Q96FV9

Purity: Affinity purification

For research use only.

Background:

HPR1 is part of the TREX (transcription/export) complex, which includes TEX1 (MIM 606929), THO2 (MIM 300395), ALY (MIM 604171), and UAP56 (MIM 142560).

To place an order, please [Click HERE](#).