

MED6

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MED6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARC33; NY-REN-28;

To place an order, please [Click HERE](#).

Catalog #: A8181

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10001

Isotype: IgG

Swiss Prot: O75586

Purity: Affinity purification

For research use only.