

SLC9A6

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 74kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC9A6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MRSA; NHE6;

Catalog #: A8187

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10479

Isotype: IgG

Swiss Prot: Q92581

Purity: Affinity purification

For research use only.

Background:

This gene encodes a sodium-hydrogen exchanger that is a member of the solute carrier family 9. The encoded protein localizes to early and recycling endosomes and may be involved in regulating endosomal pH and volume. Defects in this gene are associated with X-linked syndromic mental retardation, Christianson type. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).