

IFI44

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 50kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human IFI44

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

p44; TLDC5; MTAP44;

To place an order, please [Click HERE](#).

Catalog #: A8188

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10561

Isotype: IgG

Swiss Prot: Q8BV66

Purity: Affinity purification

For research use only.