

STK38

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STK38

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NDR; NDR1;

Catalog #: A8191

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11329

Isotype: IgG

Swiss Prot: Q15208

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the AGC serine/threonine kinase family of proteins. The kinase activity of this protein is regulated by autophosphorylation and phosphorylation by other upstream kinases. This protein has been shown to function in the cell cycle and apoptosis. This protein has also been found to regulate the protein stability and transcriptional activity of the MYC oncogene. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).