

AKR7A3

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human AKR7A3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AFAR2;

Catalog #: A8194

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22977

Isotype: IgG

Swiss Prot: O95154

Purity: Affinity purification

For research use only.

Background:

Aldo-keto reductases, such as AKR7A3, are involved in the detoxification of aldehydes and ketones.

To place an order, please [Click HERE](#).