

DAPK2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human DAPK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DRP1; DRP-1;

Catalog #: A8199

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23604

Isotype: IgG

Swiss Prot: Q9UIK4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that belongs to the serine/threonine protein kinase family. This protein contains a N-terminal protein kinase domain followed by a conserved calmodulin-binding domain with significant similarity to that of death-associated protein kinase 1 (DAPK1), a positive regulator of programmed cell death. Overexpression of this gene was shown to induce cell apoptosis. It uses multiple polyadenylation sites.

To place an order, please [Click HERE](#).