

IL36A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 18kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL36A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FIL1; FIL1E; IL1F6; IL-1F6; IL1(EPSILON); FIL1(EPSILON);

To place an order, please [Click HERE](#).

Catalog #: A8207

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27179

Isotype: IgG

Swiss Prot: Q9UHA7

Purity: Affinity purification

For research use only.