

SLC27A6

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 70kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC27A6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FATP6; ACSVL2; FACVL2; VLCS-H1;

Catalog #: A8208

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 28965

Isotype: IgG

Swiss Prot: Q9Y2P4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the fatty acid transport protein family (FATP). FATPs are involved in the uptake of long-chain fatty acids and have unique expression patterns. Alternatively spliced transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).