

HOOK2

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 83kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HOOK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HK2;

Catalog #: A8210

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29911

Isotype: IgG

Swiss Prot: Q96ED9

Purity: Affinity purification

For research use only.

Background:

Hook proteins are cytosolic coiled-coil proteins that contain conserved N-terminal domains, which attach to microtubules, and more divergent C-terminal domains, which mediate binding to organelles. The Drosophila Hook protein is a component of the endocytic compartment.

To place an order, please [Click HERE](#).