

COPS7A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 30kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human COPS7A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CSN7; CSN7A; SGN7a;

Catalog #: A8212

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 50813

Isotype: IgG

Swiss Prot: Q9UBW8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a component of the COP9 signalosome, an evolutionarily conserved multi-subunit protease that regulates the activity of the ubiquitin conjugation pathway. Alternatively spliced transcript variants that encode the same protein have been described.

To place an order, please [Click HERE](#).