

DYM

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 76kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DYM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DMC; SMC;

Catalog #: A8222

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54808

Isotype: IgG

Swiss Prot: Q7RTS9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which is necessary for normal skeletal development and brain function. Mutations in this gene are associated with two types of recessive osteochondrodysplasia, Dyggve-Melchior-Clausen (DMC) dysplasia and Smith-McCort (SMC) dysplasia, which involve both skeletal defects and mental retardation.

To place an order, please [Click HERE](#).