

MS4A12

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MS4A12

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Ms4a10;

Catalog #: A8223

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54860

Isotype: IgG

Swiss Prot: Q9NXJ0

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cell surface protein found primarily in the apical membrane of colonocytes. Silencing of this gene in colon cancer cells inhibits the proliferation, cell motility, and chemotactic invasion of cells. This gene is part of a cluster of similar genes found on chromosome 11. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).