

SEMA4C

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 93kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SEMA4C

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SEMAF; SEMAI; SEMACL1; M-SEMA-F;

To place an order, please [Click HERE](#).

Catalog #: A8225

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54910

Isotype: IgG

Swiss Prot: Q9C0C4

Purity: Affinity purification

For research use only.