

ATAD3A

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 71kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ATAD3A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8230

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55210

Isotype: IgG

Swiss Prot: Q9NVI7

Purity: Affinity purification

For research use only.

Background:

ATAD3A and ATAD3B (MIM 612317) are mitochondrial membrane proteins that contribute to the stabilization of large mitochondrial DNA (mtDNA)-protein complexes called nucleoids (He et al., 2007)

To place an order, please [Click HERE](#).