

PARL

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 42kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PARL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PSARL; PSARL1; RHBDS1; PRO2207; PSENIP2;

Catalog #: A8232

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55486

Isotype: IgG

Swiss Prot: Q9H300

Purity: Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial integral membrane protein. Following proteolytic processing of this protein, a small peptide (P-beta) is formed and translocated to the nucleus. This gene may be involved in signal transduction via regulated intramembrane proteolysis of membrane-tethered precursor proteins. Variation in this gene has been associated with increased risk for type 2 diabetes. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).