

NDUFA12

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NDUFA12

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B17.2; DAP13;

Catalog #: A8237

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55967

Isotype: IgG

Swiss Prot: Q9UI09

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which is part of mitochondrial complex 1, part of the oxidative phosphorylation system in mitochondria. Complex 1 transfers electrons to ubiquinone from NADH which establishes a proton gradient for the generation of ATP. Mutations in this gene are associated with Leigh syndrome due to mitochondrial complex 1 deficiency. Pseudogenes of this gene are located on chromosomes 5 and 13. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).