

## HS1BP3

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:100

**Calculated MW:** 43kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human HS1BP3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ETM2; HS1-BP3;

**Catalog #:** A8250

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 64342

**Isotype:** IgG

**Swiss Prot:** Q53T59

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene shares similarity with mouse Hs1bp3, an Hcls1/Hs1-interacting protein that may be involved in lymphocyte activation.

*To place an order, please [Click HERE](#).*