

CDCP1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 93kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CDCP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CD318; TRASK; SIMA135;

Catalog #: A8254

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64866

Isotype: IgG

Swiss Prot: Q9H5V8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a transmembrane protein which contains three extracellular CUB domains and acts as a substrate for Src family kinases. The protein plays a role in the tyrosine phosphorylation-dependent regulation of cellular events that are involved in tumor invasion and metastasis. Alternative splicing results in multiple transcript variants of this gene.

To place an order, please [Click HERE](#).