

RNF219

Reactivity: Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 81kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human RNF219

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C13orf7;

To place an order, please [Click HERE](#).

Catalog #: A8259

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 79596

Isotype: IgG

Swiss Prot: Q5W0B1

Purity: Affinity purification

For research use only.