

ADAM33

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 88kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ADAM33

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C20orf153; DJ964F7.1;

Catalog #: A8265

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80332

Isotype: IgG

Swiss Prot: Q9BZ11

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ADAM (a disintegrin and metalloprotease domain) family. Members of this family are membrane-anchored proteins structurally related to snake venom disintegrins, and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. This protein is a type I transmembrane protein implicated in asthma and bronchial hyperresponsiveness. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).