

ATAD3B

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 73kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ATAD3B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TOB3; AAA-TOB3;

Catalog #: A8268

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83858

Isotype: IgG

Swiss Prot: Q5T9A4

Purity: Affinity purification

For research use only.

Background:

ATAD3A (MIM 612316) and ATAD3B are mitochondrial membrane proteins that contribute to the stabilization of large mitochondrial DNA (mtDNA)-protein complexes called nucleoids (He et al., 2007

To place an order, please [Click HERE](#).