

CD200R1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CD200R1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OX2R; MOX2R; CD200R; HCRTR2;

Catalog #: A8278

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 131450

Isotype: IgG

Swiss Prot: Q8TD46

Purity: Affinity purification

For research use only.

Background:

This gene encodes a receptor for the OX-2 membrane glycoprotein. Both the receptor and substrate are cell surface glycoproteins containing two immunoglobulin-like domains. This receptor is restricted to the surfaces of myeloid lineage cells and the receptor-substrate interaction may function as a myeloid downregulatory signal. Mouse studies of a related gene suggest that this interaction may control myeloid function in a tissue-specific manner. Alternative splicing of this gene results in multiple transcript variants.

To place an order, please [Click HERE](#).