

DEDD2

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:36kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DEDD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FLAME-3;

Catalog #:A8281

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:162989

Isotype:IgG

Swiss Prot:Q8WXF8

Purity:Affinity purification

For research use only.

Background:

This gene encodes a nuclear-localized protein containing a death effector domain (DED). The encoded protein may regulate the trafficking of caspases and other proteins into the nucleus during death receptor-induced apoptosis. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).