

CEP57L1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CEP57L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

cep57R; C6orf182; bA487F23.2;

To place an order, please [Click HERE](#).

Catalog #: A8290

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 285753

Isotype: IgG

Swiss Prot: Q8IYX8

Purity: Affinity purification

For research use only.