

TICAM2

Reactivity: Human Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TICAM2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TIRP; TRAM; TIRAP3; MyD88-4; TICAM-2;

Background:

TIRP is a Toll/interleukin-1 receptor (IL1R; MIM 147810) (TIR) domain-containing adaptor protein involved in Toll receptor signaling (see TLR4; MIM 603030).

To place an order, please [Click HERE](#).

Catalog #: A8292

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 353376

Isotype: IgG

Swiss Prot: Q86XR7

Purity: Affinity purification

For research use only.