

OPRK1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human OPRK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KOR; OPRK; KOR-1; K-OR-1;

To place an order, please [Click HERE](#).

Catalog #: A8306

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4986

Isotype: IgG

Swiss Prot: P41145

Purity: Affinity purification

For research use only.