

SCAND1

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 19kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SCAND1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAZ1; SDP1;

Catalog #: A8310

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51282

Isotype: IgG

Swiss Prot: P57086

Purity: Affinity purification

For research use only.

Background:

This gene encodes a SCAN box domain-containing protein. The SCAN domain is a highly conserved, leucine-rich motif of approximately 60 aa originally found within a subfamily of zinc finger proteins. This gene belongs to a family of genes that encode an isolated SCAN domain, but no zinc finger motif. This protein binds to and may regulate the function of the transcription factor myeloid zinc finger 1B. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).