

FZD9

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 64kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FZD9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FZD3; CD349;

Catalog #: A8311

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8326

Isotype: IgG

Swiss Prot: O00144

Purity: Affinity purification

For research use only.

Background:

Members of the 'frizzled' gene family encode 7-transmembrane domain proteins that are receptors for Wnt signaling proteins. The FZD9 gene is located within the Williams syndrome common deletion region of chromosome 7, and heterozygous deletion of the FZD9 gene may contribute to the Williams syndrome phenotype. FZD9 is expressed predominantly in brain, testis, eye, skeletal muscle, and kidney.

To place an order, please [Click HERE](#).