

WNT5B

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human WNT5B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A8313

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81029

Isotype: IgG

Swiss Prot: Q9H1J7

Purity: Affinity purification

For research use only.

Background:

The WNT gene family consists of structurally related genes which encode secreted signaling proteins. These proteins have been implicated in oncogenesis and in several developmental processes, including regulation of cell fate and patterning during embryogenesis. This gene is a member of the WNT gene family. It encodes a protein which shows 94% and 80% amino acid identity to the mouse Wnt5b protein and the human WNT5A protein, respectively. Alternative splicing of this gene generates 2 transcript variants.

To place an order, please [Click HERE](#).